

## PERFORMANCE, RELIABILITY AND SECURITY



RELIABLE ENGINE AND TRANSMISSION WITH AUTOMATIC BRAKE SYSTEM



59"



27"



# DISCOVERY ORTHOQUAD .COM

Move in style with the new Orthoquad Discovery, a versatile scooter designed for everyone! Stand out with its sporty and modern look and take advantage of the latest driving technologies.

It is easy to drive and control with LCD display showing speed, battery level and distance traveled. The Discovery is also equipped with a FM - USB radio and also has a nice storage with large compartments in the steering column and under the seat.

It is equipped with a high quality suspension for a comfortable and smooth ride on any type of surface. It is the safest in its class with its fully automatic electromagnetic braking system and its manual disc brake. The Discovery has a high visibility with its complete LED lighting package for optimum security.

It's powered by a rugged, high-performance 950W engine that will never let you down and get the best warranty. The Discovery is definitely a device you can really rely on.



STORAGE UNDER THE SEAT



ALSO AVAILABLE IN SILVER

### SPECIFICATIONS

Maximum forward speed	21 km/h	Motor	950W – 24V DC
Weight capacity	150kg (330 lb)	Wheels	3-10
Turning radius	60"	Battery	12V-75h
Suspension	Front and rear	Battery charger	24V 5-amp
Range per charge	Up to 40-50km*	Braking system	Electro-magnetic fully automatic + Manual emergency brake

\*Varies with user weight, terrain type, battery charge, battery and tire condition

## SOLID CONSTRUCTION COMBINED WITH RELIABLE COMPONENTS