



**Discover a NEW WAY OF DRIVING
and enjoy your new freedom
with the Orthoquad ES400 !**



Orthoquad is proud to introduce his ES400 scooter. This all new entry-level model was designed to obtain an unbeatable value based on equipment / price ratio.

It has a long list of luxury standard equipment and features as a full lighting package including powerful LED head light, a fully adjustable deluxe orthopedic seat, a front basket, two rear view mirrors and front and rear heavy duty bumpers to insure maximum protection. It is also equipped with a delta tiller providing a better maneuverability and a safe driving experience.

Easy to drive, it is powered by a reliable 400W engine coupled with two 40 amp batteries. Each use will deliver maximum range.

**ALWAYS GET MORE
WITH THE ES400 !**



SPECIFICATIONS

Length	51"	Range per charge	Up to 30-40km*
Width	25"	Motor	400W – 24V DC
Maximum forward speed	12 km/h	Wheel size	10" pneumatic
Weight capacity	140kg (310 lbs)	Battery	12V-40Ah
Turning radius	67"	Battery charger	24V 5-amp
Suspension	Front	Braking system	Electro-magnetic fully automatic

*Varies with user weight, terrain type, battery charge, battery and tire condition

