



**Discover a NEW WAY OF DRIVING
and enjoy your new freedom
with the Orthoquad ES500 !**



Heavy Duty motor and gearbox

AVAILABLE :



4 wheel full suspension

The all-new **ES500** from Orthoquad will give you comfort, safety, reliability and increased efficiency. This new model is equipped with a sturdy 450W high quality engine that ensures peace of mind and exceptional reliability. Paired with two 40 amp batteries, it offers an exemplary autonomy with each use. Its 4-wheel dynamic suspension delivers optimum performance, a smooth and comfortable ride and superior stability on all types of terrain.

With the **ES500** you will find an ergonomic design and intuitive operation for safe movement and easy use. In order to provide you with a maximum experience, the **ES500** comes with an impressive list of standard features such as a high-brightness LED headlamp, a fully adjustable rotatable luxury orthopedic seat, pneumatic wheels, closed front basket, two mirrors and sturdy front and rear bumpers for maximum protection. It is also equipped with a delta-type steering wheel for better maneuverability and safe driving.

The **ES500** is easy to drive and designed to meet the essential needs of new or experienced users.

**The perfect combination of
comfort, safety and reliability!**

SPECIFICATIONS

Length	51"	Range per charge	Up to 30-40km*
Width	25"	Motor	450W – 24V DC
Maximum forward speed	12 km/h	Wheel size	10" pneumatic
Weight capacity	140kg (310 lbs)	Battery	12V-40Ah
Scooter weight W/O batt	75kg (165 lbs)	Battery charger	24V 5-amp
Suspension	Front	Braking system	Electro-magnetic fully automatic

*Varies with user weight, terrain type, battery charge, battery and tire condition

**AGM or gel type required

