

GS1400

OrthoQuad is proud to introduce it's new GS1400, which has nothing to compare with. The GS1400 is nothing else than one of the most powerfull and comfortable mobility scooter. With its 14 inch wheels, nothing can stop it. Its incredible four wheels suspension is built to ride smoothly even on rough tracks. If you like to ride around fields or lakes, the GS1400 is definitely the mobility scooter you deserve.

The GS1400 is alone in its class. A road test will convince you.



Heavy duty model, max load : 450 pounds


A canopy and deluxe rear basket are includ with the GS1400



RED



GREEN



Specifications

Length	62"	Charger	120V - 5A
Width	28"	Batteries	2 x 12V - 75Ah (group 24)
Max. Speed	16 km/h	Wheels	14 "
Range (Up to)	30 - 40 km	Turning radius	1.65 m - 5,4'
Weight capacity	205 kg - 450 lbs	Braking distance	1.8 m - 6'
Suspension	full active front and rear	Drive Train	1400W - 4 Pole
Max. climbing gradient	15 degrees		