



ORTHOQUAD

New Generation .COM

MONARCH

The new generation Monarch wheelchair establishes a new bench mark in terms of design and performance. Built with a high level active suspension, it is designed for inside and outside use. The Monarch's drive train is equipped with two motors for a maximum of torque and power. With its effortless ride and handling, the Monarch makes your mobility accessible and satisfying.



Specifications

Weight capacity	280 lbs - 125kg	Rear wheels	6"
Speed	6.5 km/h	Drive train	2 motors, mid wheel drive
Width	24.5"	Braking system	Electronic, regenerative disc brakes
Length	44"	Electronic	PG VR2 40A
Drive wheels	10"	Charger	Off-board 5A
Front wheels	5"	Range	Up to 25 km
		Battery	12V - 33Ah - U1 x 2

