

RS1200

Here is the new RS1200 from ORTHOQUAD. It redefines the standards of comfort and handling. The conception of this revolutionary scooter has been designed by high level engineering. The RS1200 is a high quality mobility scooter with amazing performances. It has a powerful 1200 watt engine coupled to a heavy duty transmission. Acceleration and automatic braking system are perfectly working together.

With its big 14 inch wheel and an improved 4 wheel independent suspension, the RS1200 will give you an incredibly smooth ride.

Definitely, the RS1200 will satisfy anyone



Silver



Dark green



Blue



Specifications:
 Length: 58"
 Width: 27"
 Max. Speed: 15 km/h
 Range(approx.): Jusqu'à 50 km

Max. Load: 180kg - 400lbs
 Max. climbing gradient: 15°
 Engine: 1200W
 Suspension - Front and rear
 Turning radius: 1.65m - 5,4'

Charger: 120V - 5Ah
 Batteries: 2 x 12V - 75A (Included)
 Wheels: 14"
 Braking distance: 1.8m - 6'
 Weight: 122kg - 268lbs

