




RS1400
19 KM/H EDITION

Now with heavy duty 8A charger



Orthoquad is proud to introduce its new generation full size scoter, the RS1400. Designed and developed for heavy duty performance, it also includes all the luxury features that have come to be expected from the market. A long list of standard equipments as a full lighting package, deluxe high back seat and delta tiller for easy throttle operation with fingers provide a safe and practical mobility experience. The RS1400 delivers a premium performance with a 1400 watt four pole motor for exceptional power, torque and efficiency.

Excellent for outdoor performance with its 14" all-terrain tires, amazing ground clearance and a 4 wheel full suspension, the RS1400 will give you an optimal level of comfort on any surface. With the RS1400, you will find speed, comfort, range and reliability.

The RS1400 is definitely your best choice!



DS1400
Dual seat edition

SPECIFICATIONS

| | |
|--------------------------------|----------------------------|
| Length | 63" |
| Width | 27" |
| Max. Speed | 19 km/h |
| Range (Up to) | 30 - 40 km |
| Weight capacity | 180 kg - 400 lbs |
| Suspension | full active front and rear |
| Transport weight W-O batteries | 85 kg - 187 lbs |

| | |
|------------------------|---------------------------|
| Charger | 120V - 8A (New) |
| Batteries | 2 x 12V - 75Ah (group 24) |
| Wheels | 14" |
| Turning radius | 1.65 m - 5,4' |
| Braking distance | 1.8 m - 6' |
| Drive Train | 1400W - 4 Pole |
| Max. climbing gradient | 15 degrees |