



Discover a **NEW WAY OF DRIVING**  
and enjoy your new freedom  
with the Orthoquad **RZ1000!**

DIGITAL DISPLAY



Bluetooth  
speaker



USB charging  
ports



FM Radio



All new Orthoquad RZ1000 establishes a new bench mark for mid-size scooters. Powered by the reliable and powerful high performance 950 watt engine and coupled to 2x 75Ah batteries, it delivers a 15km/h top speed and an impressive range. The RZ1000 offers a high level of comfort even on rugged terrain with its full 4 wheel suspension.

It also has a long list of luxury standard features: complete LED lighting package for maximum brightness performance and minimum maintenance, new fully adjustable deluxe high-back seat with headrest and sliders, lockable front basket, bottle holder, two rear view mirrors, front and rear heavy duty bumpers to insure maximum security. The RZ1000 also has an FM radio with bluetooth speaker and its USB ports allows to charge your small electronic devices such as your mobile phone.

The new redesigned tiller has a gas sprung activator providing infinite tilt steering adjustment for ergonomic user positioning. Improved delta tiller throttle control selection can be activated by fingers or thumb of both hands. Get the latest driving assistance technology with the new easy to read control panel with integrated digital display.

**The RZ1000 will take you  
almost anywhere with style  
and safety.**



Metallic Blue

Apple Candy

#### SPECIFICATIONS

Length	56"
Width	26"
Maximum forward speed	15 km/h
Weight capacity	150kg (330 lbs)
Turning radius	RZ1000 : 67"
Suspension	Front and rear

Range per charge	Up to 40km*
Motor	950W – 24V DC
Wheel size	110/55-8
Scooter weight W/O batteries	110kg (242lbs)
Battery	12V-75Ah**
Battery charger	24V 5-amp
Braking system	Electro-magnetic fully automatic

\*Varies with user weight, terrain type, battery charge, battery and tire condition

\*\*AGM or gel cell type required

