



**Discover a NEW WAY OF DRIVING
and enjoy your new freedom
with the Orthoquad RZ1200 !**

Orthoquad is proud to introduce its new full size scooter, the **RZ1200**. Designed and developed for heavy duty performances, it also includes all the luxury features that have come to be expected from the market. A long list of standard equipments as a full LED lighting package for maximum brightness, performance and minimum maintenance, deluxe high back seat with headrest and sliders, lockable front basket, bottle holder, two rear view mirrors, front and rear heavy duty bumpers to insure maximum security, and now with two USB charging ports for small electronic devices as cell phones.

The new redesigned tiller has a gas sprung activator providing infinite tilt steering adjustment for ergonomic user positioning. Its easy operation provides a safe and practical mobility experience. Get the latest driving assistance technology with the new easy to read control panel with integrated digital display. The **RZ1200** delivers a premium performance with our legendary heavy duty 950 watt four pole motor for exceptional efficiency.

The **RZ1200** is excellent for outdoor performances with its 14" all-terrain tires, amazing ground clearance and 4 wheel full suspension. It will give you an optimal level of comfort on any surface.

**With the RZ1200, you will find speed,
comfort, range and reliability.
The RZ1200 is definitely your best choice!**



SPECIFICATIONS

Length	63"	Range per charge	Up to 40km*
Width	27"	Motor	950W – 24V DC
Maximum forward speed	16 km/h	Wheel size	4.00-6 / 14"
Weight capacity	180 kg (400 lbs)	Battery	12V-75Ah** gr.24
Scooter weight W/O batt	110 kg (242 lbs)	Battery charger	24V-5Amp
Suspension	Front and Rear	Braking system	Electro-magnetic fully automatic

*Varies with user weight, terrain type, battery charge, battery and tire condition
**AGM or gel type required

