

RZ1500-R ORTHOQUAD .COM

Discover a **NEW WAY OF DRIVING**
and enjoy your freedom with
the Orthoquad **RZ1500-R!**



Bluetooth speaker



USB charging ports



FM Radio



Heated seat



AVAILABLE:



RED



BLACK



CHARCOAL
GREY



SPECIFICATIONS

Length	67"
Width	31"
Maximum forward speed	21 km/h
Weight capacity	160 kg (350 lbs)
Scooter weight W/O batt	130 kg (286 lbs)
Suspension	Front and Rear

Range per charge	Up to 40-50km*
Motor	1400W – 24V DC
Front wheel size	14" pneumatic (14x4.00-8)
Rear wheels size	16" pneumatic (130/60-10)
Battery	12V-100Ah
Battery charger	24V 8-amp
Braking system	Electro-magnetic fully automatic

The new Orthoquad **RZ1500-R** is our fastest and most luxurious mobility scooter ever made. It is the result of several years of development in order to provide an unequaled equipment/performance ratio. Powered by its reliable and powerful 1400 watt motor and paired with two 100Ah batteries, it can reach an impressive speed of 21 km/h, with great range. With its four wheel active suspension, the **RZ1500-R** offers a high level of comfort even on rough terrain.

It has a long list of standard luxury features such as full set of LED lights, a fully adjustable high-back heated seat with adjustable headrest and armrests, a lockable front and rear basket, a cup holder, two mirrors and solid front and rear bumpers ensuring maximum safety. It is also equipped with a set of front and rear mudguards for optimal protection. The **RZ1500-R** also has an FM radio with bluetooth speaker and its USB ports allows to charge your small electronic devices such as your mobile phone.

The tiller is equipped with an easily adjustable hydraulic system for ergonomic positioning of the user. It also features an upgraded delta steering wheel that allows hand or thumb throttle control. You will also find the very latest driver assistance technologies in the new easy-to-read control panel with integrated digital display showing your speed, battery status, time and even the outside temperature. The **RZ1500-R** is protected by a 2 year engine and transmission warranty.

OUR FASTEST AND MOST LUXURIOUS MOBILITY SCOOTER!



*Varies with user weight, terrain type, battery charge, battery and tire condition
**AGM or gel type required