



Discover a **NEW WAY OF DRIVING**  
and enjoy your new freedom with  
Orthoquad **RZ1500**

DIGITAL DISPLAY



Bluetooth  
speaker



USB charging  
ports



FM Radio



Electric blue



Charcoal Grey

The all new modern RZ1500 raises the bar for luxury, comfort and performance. Powered by our renowned reliable and powerful 1400W engine and coupled to two 100 amps batteries, it has a top speed of 18 km/h and an impressive range.

Totally redesigned and developed for dynamic performance in every condition, the RZ1500 delivers amazing standard features like LED lighting package for maximum brightness performance and minimum maintenance, bottle holder, front lockable basket, two rear view mirrors, 8A heavy duty charger. The RZ1500 also has an FM radio with bluetooth speaker and its USB ports allows to charge your small electronic devices such as your mobile phone.

Get more comfort with the new tiller providing infinite tilt steering adjustment for ergonomic user positioning. The RZ1500 is also equipped with a standard fully adjustable deluxe high-back orthopedic seat.

Get the latest driving assistance technology with the easy to read control panel with integrated digital display.

**Get what you deserve  
with the new RZ1500!**



## SPECIFICATIONS

Length	67"
Width	27"
Maximum forward speed	18 km/h
Weight capacity	180kg (400 lbs)
Turning radius	71"
Suspension	Front and rear

Range per charge	Up to 40-50km*
Motor	1400W – 24V DC
Wheel size	14 x 4.00-8
Scooter weight W/O batt	125kg (275lbs)
Battery	12V-100Ah**
Battery charger	24V 8-amp
Braking system	Electro-magnetic fully automatic

\*Varies with user weight, terrain type, battery charge, battery and tire condition

\*\*AGM or gel cell type required

