



**Discover a NEW WAY OF DRIVING
and enjoy your new freedom with
Orthoquad RZ1500**

DIGITAL DISPLAY



The all new modern RZ1500 raises the bar for luxury, comfort and performance. Powered by our renowned reliable and powerful 1400W engine and coupled to two 100 amps batteries, it has a top speed of 18 km/h and an impressive range.

Totally redesigned and developed for dynamic performance in every condition, the RZ1500 delivers amazing standard features like LED lighting package for maximum brightness performance and minimum maintenance, bottle holder, front lockable basket, two rear view mirrors, 8A heavy duty charger and two USB charging ports for small electronic devices as cell phones.


Get more comfort with the new tiller providing infinite tilt steering adjustment for ergonomic user positioning. The RZ1500 is also equipped with a standard fully adjustable deluxe high-back orthopedic seat.

Get the latest driving assistance technology with the easy to read control panel with integrated digital display.

**GET WHAT YOU DESERVE
WITH THE NEW RZ1500!**



Available:



SPECIFICATIONS

Length	67"	Range per charge	Up to 40-50km*
Width	27"	Motor	1400W – 24V DC
Maximum forward speed	18 km/h	Wheel size	14 x 4.00-8
Weight capacity	180kg (400 lbs)	Battery	12V-100Ah**
Turning radius	71"	Battery charger	24V 8-amp
Suspension	Front and rear	Braking system	Electro-magnetic fully automatic

*Varies with user weight, terrain type, battery charge, battery and tire condition
**AGM or gel cell type required

