



We don't believe mobility scooters should be about compromise. We believe your mobility scooter should improve your life. That is why GNX introduced the 4 season Orthoquad Summit.

The summit keeps you protected for a comfortable ride in all weather's conditions. It is definitely the most advanced and unique mobility scooter on the market!

Featuring reverse camera, parking sensors, FM Radio, USB charging ports and internal heater, it really provides comfort's level of a car. The soft suspension and 100Ah long range batteries will keep you going for hours. And its full digital dashboard makes driving simple!



DIGITAL DASHBOARD



REVERSING CAMERA



LUXURY CAPTAIN SEAT



PLUG IN CHARGING SYSTEM



INTERIOR HEATER



WINDSHIELD WASHER



LED HEAD LIGHTS



REMOTE CENTRAL LOCKING



REAR LED LIGHTS



REAR PARKING SENSORS



ADVANCED SUSPENSION



Smooth running electric motor



SPORT ALLOY WHEELS



e-mobility

SPECIFICATIONS

Length	66"
Width	30" (Without mirrors)
Height	66"
Maximum forward speed	15 km/h
Weight capacity	160kg (350 lbs)
Scooter Weight	245kg (540lbs) With batteries

Suspension	Front and rear
Range per charge	Up to 40km*
Motor	1400W - 24V DC
Wheel size	14" pneumatic (14x4.00-8)
Battery	12V-100Ah**
Battery charger	24V 8-amp
Braking system	Electro-magnetic fully automatic

*Varies with user weight, terrain type, battery charge, battery and tire condition
 **AGM or gel cell type required

